

Toilet Support Rail

The **3 in 1 Rail** is designed for safety and comfort and benefits a range of users from the average build to the larger build persons. This Rail enables the user to have support close to the body, where it is most needed. This model has a weight loading of 110Kg. The user is able to lower safely onto the toilet and pushing up from the seat is also much easier. The Rail is safer than over frames as it is rigidly fixed to the toilet pan, and also does not need to be removed for other users. The handles can also be folded down in order to facilitate a wheelchair transfer. Additionally, either one or both handles can be removed to suit some requirements.



Standard Mode

Splayed Mode

Fold Down Mode

RECOMMENDED FOR:

Multipurpose use

SPECIFICATIONS:

Two optional settings: 510mm or 640mm between the handles, and fold down handles
Height of handles above seat 260mm approx
The back plate internal measurements 370mm



Easy to activate emergency buttons can be fitted to Rails for safety and convenience.